

ALTIKA

SCENAR-therapy (Self-Controlled Energy Neuro Adaptive Regulation)



SCENAR-therapy (Self-Controlled Energy Neuro Adaptive Regulation)

A series of unique physiotherapeutic equipment having biological feedback for treatment and prevention of a large number of diseases by activating the human neuroregulating system. SCENAR-therapy uses the newest achievements in dynamic medicine, micro-electronics, microprocessor engineering.

SCENAR-therapy

Efficient drug-free Self-Controlled Energy Neuro Adaptive Regulation of an organism

"A Combination of Western and Eastern Medicine is similar to divine Love being an ideal of perfection"

Numie Yamashita(I.Nakatani, N. Yamashita, "Riodaraki Acupuncture", 1977)

The SCENAR-therapy is carried out with the help of SCENAR equipment. It is an extra small electronics with the verified methodical support, friendly interface and weight of only 300 grams.

SCENAR equipment differs from conventional electrostimulators due to availability of biological feedback: the shaping of electropulses is controlled by the electroskin resistance of the patient changing under effect of these pulses.

SCENAR equipment is a polyclinic, first aid, and mini-hospital in your pocket

SCENAR-97.0. "The Home physician". Modern version of the equipment with mode display and map of response dynamic in an individually - dose mode. Spot frequencies: 14, 60, 120, 320 Hz, intermittent duty 3:1; a sweep frequency mode.

SCENAR-97.1. Modern version of professional equipment operating in continuous and modulated modes. The digital control allows to change the frequency from 15 to 350 Hz; to provide continuous mode of operation with a ratio of 1:1; 2:1; 3:1; 4:1; 5:1; to influence in a sweep frequency mode of 30 to 120 Hz. The mode display is made using LED and sound.

SCENAR-97.2. The same performances as SCENAR-97.1. Besides it has an individually - dose mode (IDM) and

indication of response parameters, that allows to optimise the therapy.

SCENAR-97.3. Has the same functions as SCENAR-97.1. Liquid-crystal mode display. The "intense" therapy mode with the option of adjustment of distance between stimulant pacemakers (I = 2...8).

SCENAR-97.4. Has the same functions as SCENAR-97.3. Liquid-crystal mode display and IDM. The response parameters (on the LCD) allow to optimise therapy.

SCENAR-97.7. The computer version including the SCENAR-97 equipment, module of computer interface unit and software packages: information (IS), real time (RTS), a posteriori processing (APPS). IS contains methodical software with the instructions demonstrating operation in real time when treating specific diseases. The opportunity of computer storage of specially developed methods. RTS allows to input, to process and to store the SCENAR impact during treatment. APPS ensures presentation of the stored information in convenient form.

All the equipment is patented and certified in Gosstandart of Russia, certificate ? 02754386.

Operation principles

SCENAR - effect activates a large part of nerve fibres, including thin C-fibres. It creates an efficient dose of regulator peptides and other humoral compounds that are capable to manage with a broad spectrum of diseases that are caused not only by function disorders, but also processes with expressed tissue, cellular and subcellular damages. Thus endophrynum removes pain; vasoactive intestinal peptide, somatostatynum and bombezynum normalise the gastrointestinal secretion; vasopressynum and adrenocorticotropynum actuates learnability and memory; luliberynum influences sexual activity; somatostatynum depresses temperature; thyroliberinum improves respiratory function, etc.

Effects of SCENAR-therapy

- /// Normalisation of the infringed functions,
- /// Acceleration and reduction of pronouncement of pathological processes
- /// Normalisation of exchange processes
- /// Regulation of vascular tone
- /// Considerable improving of general condition with increase of adaptive resources of an organism

- /// Improving the psycho-emotional condition
- /// Sedative, anti-edemic, anti-allergen, antiphlogistic effects

The first procedures give relief. The procedures are carried out on a daily basis or 2-3 times a week depending on disease. The mean course of treatment consists of 8 to 10 sessions. When curing chronic diseases the courses can be repeated.

Advantages of SCENAR-therapy

- /// Versatility. Can be used for treating different diseases in different conditions
- /// Variability. Enables to change treatment parameters
- /// Biological controllability. The parameters of application are adjusted by organism itself in each specific instant
- /// Physiological adaptability. Ensures a selective absorption of applied energy in the amount, needed for correcting the damaged functions
- /// Non-invasive effect. The technique is simple and harmless
- /// Can be applied in any age
- /// High efficiency

The indications and contra-indications

The application of SCENAR-therapy is indicated at any stage of diseases:

- /// Nervous system - (Neuralgia, neuritis - impairment of cerebral circulation and its consequences, vegetative neuralgia, neurosis, trauma and inflammation of central nervous system, different spine diseases with secondary disorders of nervous system, etc.)
- /// Osseous-muscular system - (myositis, arthritis, arthrosis, bruises of soft tissues, fractures)
- /// Respiratory tract - (disorders of the upper respiratory tract, tracheitis, bronchitis, pleurisy, pneumonia, bronchial asthma)
- /// Cardiovascular system - (ischemic heart disease, hypertensive disease, cardiac arrhythmia, lower extremities circulation disease, varicose diseases, trophic ulcers, etc.)
- /// Urogenital system - (cystitis, pyelonephritis, adnexitis, prostatitis, renal colic, disordered menstrual cycles, impotence, frigidity, etc.)
- /// Dental disorders - (parodontosis, parodontitis, pulpitis, etc.)

- /// and other pathologies.

Contra-indications

- /// Individual intolerance
- /// Availability of the artificial rhythm pacer
- /// Heavy mental diseases
- /// Acute infectious diseases with obscure diagnosis.

Who can use SCENAR

Today users of the equipment in Russia are physicians and medical staff of virtually all the licensed specialities, including neurologists, therapists, surgeons of different specialities, traumatologists, gynaecologists, accoucheurs, dermatologists, otolaryngologists, stomatologists, sporting doctors, physiotherapists, experts in resuscitation and so on.

The doctors and experts are capable to render the qualified assistance to patients in any situations.

The special version of the equipment is accessible for persons without medical education to use at home and office, as well as at water transport, in automobiles, etc.

Whether the training is required?

Doubtlessly, for the professionals meeting serious challenge, the training is necessary. It is carried out by the leading experts on SCENAR-therapy. At the same time, developed methodical maintenance (the operation manual contains techniques of treatment of many diseases) and absence of contra-indications to SCENAR-therapy, allow to use it in simple situations, after initial training at specialised centres. The users are not obliged to know the topography of acupuncture points or to conduct their search. 80% of success is due to effect on the projection of the pathological centre and its symmetric zones. The time of procedures makes about 10 minutes in simple situations.

Whether the price is efficient?

The economies of application of SCENAR are invaluable. Its use helps not only to essentially reduce the consumption of pharmaceuticals, to increase the protective properties of an organism, to moderate the number of visits to the doctor, to be relief pain, but also to help patients in situations, when other techniques do not help.

For more information, contact:

Lars Harrison, CEO, Altika

6325-9 Falls of Neuse Rd., Suite 193, Raleigh, NC 27615-6809, USA

Phone: 919-878-8166, Cellular: 919-523-1770, Efax: 651-324-3440

lharrison@altika.com

Program code: 404-SCENAR